

LABEYE

Designed specifically for laboratory use, the LabEye is a ready-to-implement multi-camera system performing visual process monitoring, chronicling research in HD clarity and providing automation troubleshooting. Engineered by Two Square Science this collaborative effort between Distek, Two Square Science and Industry leaders offers a welcome and needed addition to any laboratory that looks to streamline pharmaceutical testing in both research and development and quality environments.

- **LAB READY CAMERAS**

Capturing images in true HD, the LabEye cameras record 30fps at 1080p resolution.

- **DIRECT TO DIGITAL, HD RESOLUTION**

Take snapshots at a high frame rate with HD resolution and watch live or save video for later playback.

- **SCHEDULING**

Timed snapshot and video recorded events can be scheduled and created in the control software then sent to the command center which will start processing the schedule immediately.

- **STANDARD MOUNTING SYSTEM**

Equipped with a standard mount system for popular dissolution systems and other analytical devices.

- **SCALABLE**

Scalable from 1 - 8 cameras.

- **FILE VALIDATION**

All files are stored along with a unique file name and key computed from the data in the file. A file moved from the Command Center can be re-uploaded and validated to ensure file edits have not occurred.

- **RAPID ROI**

Low cost of ownership and implementation with a rapid return on investment. Free up high value personnel.



STANDARD SYSTEM PACKAGE

INCLUDED	
LabEye Cameras	<p>The LabEye Cameras are true HD and record 30fps at 1080p resolution.</p> <ul style="list-style-type: none"> • 1 - 8 Cameras • Protective Lab Housings • Connection Cables
LabEye Brackets	Based Upon Customer Requirements
LabEye Command Center	<p>The LabEye Command Center connects to a Laptop or PC via a LAN port and is equipped with up to 8 USB ports depending on the camera configuration.</p> <ul style="list-style-type: none"> • 2Tb (Terabyte) Hard Drive • Power Adapter - 110V AC Power Connection <p>(Unless otherwise specified. Please specify country of destination.)</p>
LabEye Control Laptop	Windows 7 Professional



Distek symphony 7100
Bathless Dissolution
LabEye Bracket Configuration



Distek Model 2500
Water Bath Dissolution
LabEye Bracket Configuration